


LesionQuant

NeuroQuant

Quick Start Guide

Installed System



**For LesionQuant users on an
Installed System.**

Follow this checklist to
correctly set up and process
3D T1 and FLAIR MR images.



CORTECHS Labs

Better insight through quantitative imaging.

Start

Getting Started

- ❑ Download recommended scanner settings for both NeuroQuant and LesionQuant for GE, Philips, Siemens, Toshiba (1.5 T and 3.0 T) and Hitachi (1.2 T, 1.5 T and 3.0T)
- ❑ Ensure that NeuroQuant can be run on your operating system. Supported systems include Mac OS X, Microsoft Windows and Ubuntu Linux.
 - ❑ See Hardware Requirements for detailed information. www.cortechslabs.com/resources/installed-system

Logging In

- ❑ Using the default User ID [admin] and Password [ctxapp], log into NeuroQuant: <https://hostname>
 - ❑ Replace 'hostname' with the IP address of your NeuroQuant system, if logging in from a different computer.
- ❑ After logging in, you will be directed to change your password.
- ❑ Your password* must have:
 - ❑ At least eight (8) characters
 - ❑ At least one UPPER CASE letter
 - ❑ At least one special character (e.g. @,#,\$,%, etc.)
- ❑ Confirm new password by entering twice.

* Protect the privacy of your account and data by using a strong and hard to guess password.



Site Configuration/System Admin

- Configure site specific settings on the *Site Configuration* page.
 - Enter the correct administrator email and valid SMTP server (if required) to ensure system messages are received.
 - Modify the data retention period from the default (365 days), if desired.
 - Select the minimum lesion size and separation values, as needed.
- Click Submit to save changes. All changes will be lost if submit is not clicked before navigating away from the page.

NOTE: DO NOT change the Site Name or your license key will be invalidated.


Queue Monitor

Dark Background

- Checking the Dark Background box will convert all reports to a dark background.
- Even when the Dark Background box is checked, white background PDF reports are available by download in the *Queue Monitor* page.

CSV Output File

- Must be enabled by license to use.
- Clicking the CSV icon will export all data in a user account as a ZIP file.



Process

Process 3D T1 and FLAIR Scans

Using port 104 on the NeuroQuant system:

- ❑ Select the destination DICOM AETitle that corresponds to the required report.

Report to Generate	Destination DICOM AETitle
LesionQuant FLAIR Report	lesionquant
LesionQuant FLAIR Report PLUS	lesionquant_plus

- ❑ If two time points are available, send the **prior/baseline series** images to the same report type, matching the report desired (if not, skip this step):
 - ❑ T1-weighted, sagittal, non-contrast enhanced 3D MR series as DICOM images acquired using NeuroQuant recommended scanner settings
 - ❑ Matching 2D or 3D FLAIR axial or sagittal series as DICOM images
 - ❑ Ensure that PatientID and StudyInstanceUID matches in this uploaded 3D T1/FLAIR pair.
- ❑ Send the **new/current series** of images to the same report type, matching the report desired:
 - ❑ T1-weighted, sagittal, non-contrast enhanced 3D MR series as DICOM images
 - ❑ Matching 2D or 3D FLAIR axial or sagittal series as DICOM images
 - ❑ Ensure that PatientID and StudyInstanceUID matches in this uploaded 3D T1/FLAIR pair
- ❑ When the upload dialog box closes, your series will be automatically queued for processing.
 - ❑ It may take several minutes for your series to be transferred.
- ❑ Check that upload progress is displayed.
 - ❑ Processing will start automatically.
 - ❑ It could take up to 30 minutes to process, depending on computer performance specifications.

LesionQuant Complete

- ❑ Output will be sent automatically in DICOM format to the PACS host configured in the *System Administration* settings.
- ❑ Consult the NeuroQuant user manual or support@cortechslabs.com for more information.



Finish



CORTECHS Labs

For more tips and information on LesionQuant, visit our [Cortechs Labs Blog](#)